Device comprising a support base and a handle or pestle equipped with a force indicator

3 interchangeable application heads 3 different shapes and sizes





A protective cover to adapt on the base of the MyoDK

> A suitcase for storage and transport





RK INNOV SAS 25 rue des écoles, 95500 Le Thillay (France) SIRET 83428848200021 GLN 3770018797014 SRN N° FR-MF-000028533

> Contact contact@rkinnov.fr www.pcptherapy.com

RK INNOV products

Medical and paramedical devices, designed and developed by RK INNOV worldwide patented and made in France Company RK INNOV certified EN ISO 13485

RK INNOV trainings

Practice of the technique PCPTherapy Use of the MyoDK device for PCPTherapy Use of the MyoRK equipment for vacuum therapy



La certification qualité a été délivrée au titre de la ou des catégories d'actions suivantes : ACTIONS DE FORMATION



A Class I medical device, effective and easy to use



A mechanical device, specifically designed for the application of positive pressures according to the PCPTherapy technique*, by using the body force. * see PCP Brochure

It is intended to relieve the musculoskeletal pain and to prevent and treat pathologies of the locomotor apparatus.

It allows to treat the whole body in one session.

Technical novelty

- Use of the body weight instead of hand force allows to relieve upper limbs and reduces the therapist's fatigue in his daily repetitive actions
 - Indication of the level of the applied force to control the applied pressure
 - Patented application heads with multiple application points.

MDK Brochure E 231228.docx

Use performance

- Portable device, storage in a lockable case
- Simple manual use, no setting
- Comfortable use: design with a support base, equipped with a protective cover, to use body weight and thus relieve his upper limbs
- Use flexibility in different positions and application methods allowing multiple uses for personalized and adapted treatments
- 3 interchangeable application heads (TAMs)
 - Quick connection and reliable connection
 - Ergonomics for a comfortable application
 - multiple pins working simultaneously for quick work on the whole body
- Possibility to combine PCPTherapy with the MyoDK device and vacuum therapy with the MyoRK equipment (see MRK Brochure)

Performance requirements

- Applied force: up to 250 N (25 kgf)
- Session duration: 20 min

Design safety

- Mechanical device offering most of the advantages of manual techniques
- Indicator displaying the level of the applied force allowing to determine the pressure pain threshold and to control the applied pressure
 Application heads
- Different ergonomic shapes and sizes allowing to respect the different anatomical areas
- Equipped with several and wide pins allowing to both have multiple application points and a lower pressure on skin-contact surfaces and then to reduce the pain felt at the pressure
- Smooth surfaces for ease of decontamination
- Biocompatible material to avoid any biological risks

Clinical benefits

Relieve acute and chronic pain from the first sessions

- Reduction of the myofascial and tendinous tensions (general relaxation)
- Relaxation of muscles, including deep muscle
- Decrease of joint constraints
- Reduction of the pain
- (Increase of the pressure pain threshold)
- Increase of joint flexibility and mobility
- Posture improvement
- Quality of life improvement

Indications

- Back pain including low back pain and neck pain and joint pain
- Prevention, supporting and/or treatment of the
- pathologies of locomotor apparatus
 - Tendinitis,
 - Rachialgia (cervicalgia, dorsalgia, lombalgia),
 - Posture disorders
 - Fibromyalgia, ...
- Sport



Use safety

Use limited to the health professionals specialized in manual therapies after a specific training provided by RK INNOV (see the brochure on www.pcptherapy.com)







2007 : Lauréat Concours National Ministère de l'enseignement supérieur et de la recherche





Partenariat & collaboration Commissariat Energie Atomique





14 international patent families Find out more about <u>www.pcptherapy.com</u>